



# Chait

PERCUTANEOUS CECOSTOMY CATHETER




Improve quality of life for your patients.



# Chait


## PERCUTANEOUS CECOSTOMY CATHETER



STIFFENING CANNULA




CHAIT CECOSTOMY CATHETER  
Polyurethane



CHAIT ACCESS ADAPTER WITH CONNECTING TUBE

- Placement is a minimally invasive, reversible procedure.
- Facilitates antegrade colonic cleansing for patients with fecal incontinence.
- Unique, low-profile trapdoor fitting designed to lie flat against the skin.
- Enables greater patient independence to manage fecal incontinence.
- Catheter is 10.2 Fr and features three sideports.



DAWSON-MUELLER MULTIPURPOSE DRAINAGE CATHETER

### DAWSON-MUELLER MULTIPURPOSE DRAINAGE CATHETER

- Mac-Loc® locking loop helps prevent catheter displacement.
- Soft, smooth Ultrathane® material helps maintain patient comfort.
- Hydrophilic coating eases insertion.

Order Number	Reference Part Number	Tract Length cm	Extended Catheter Length cm
<b>Chait Percutaneous Cecostomy Catheter</b>			
G11632	TDCS-100	0-6	20
G36392	TDCS-100-M	3-9	24
G12685	TDCS-100-L	6-14	28
<b>Chait Percutaneous Cecostomy Set – Including Short-length Chait (TDCS-100) and Accessories Below</b>			
G12571	TDCS-100-CS		

## ACCESSORIES

### Dawson-Mueller Multipurpose Drainage Catheter

G09705 ULT8.5-38-25-P-5S-CLDM-HC

### Chait Access Adapter with Connecting Tube

G11630 CTU10.0-18-CAA

### Chait Enterostomy Suture Anchor Set

G12230 GIAS-100-CHAIT

### Amplatz Ultra-Stiff Wire Guide

G10168 THSF-35-80-AUS

### Coons Dilator

G03927 JCD8.0-38-20-COONS

- AORTIC INTERVENTION
- CRITICAL CARE
- ENDOSCOPY
- INTERVENTIONAL RADIOLOGY
- LEAD MANAGEMENT
- PERIPHERAL INTERVENTION
- SURGERY
- UROLOGY
- WOMEN'S HEALTH

### Customer Service Centers

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

