



PRODUCT INFORMATION

Product-No.: 50.02.36 – (No. non-sterile: 04.47.)
Name: Combi Cap Luer lock male Luer lock female - PP – red – s t e r i l e
Production: The product will be produced under a certified QMS according to EN ISO 13485

The following materials are used:

Product	Material
Combi Cap Colour	01.11. Purell HP 548 N 01.25. MB PE 3927 - red

Packaging: upper foil: paper
lower foil: PVC-foil (DEHP-free)

Packaging around: white dispenser carton: 12,5 x 9 x 7,2 cm - corrugated cardboard
brown carton: 37,5 x 25,8 x 46,5 cm - cardboard

Substances of animal origin: Tallow derived additives may be used in the manufacture of Purell HP 548 N. Tallow derived materials used in this product fulfill the requirements laid down in the Regulations 1069/2009/EC, and 142/2011/EC, and the "Note for Guidance EMEA/410/01, rev. 3".

REACH: The product doesn't contain any danger material following REACH

94/62/EG: Product and packaging contain in the sum less than 100 ppm of lead, cadmium, mercury and 6-valued chrome

Sterilization: The sterilization is effected by ethylene oxide.
The sterility of the products and also the efficiency of the sterilization process are tested by an accredited laboratory.

Risk assessment: The visible risks as incompleteness or dirtiness are minimized by suitable checks and measures. Therefor an endangering of the patient can be excluded in case of an appropriated use of the product.

The product doesn't contain any latex.